

Lite Wand LED
User Manual

Thompson Retractor 



Uncompromised Exposure

Utilize the bright, reusable, Lite Wand LED to add illumination to any open surgical procedure with the Thompson Retractor.



STEP 1

Attach joint to any 1/2" frame component where lighting is desired.



STEP 2

Clip the Lite Wand LED to the joint and retract to the desired position.



STEP 3

Secure the Lite Wand LED by tightening the joint.



STEP 4

Plug in the Lite Wand LED cable to the Thompson LED light source. Turn on light source and position the flexible neck of the Lite Wand LED as needed.

NOTE: Lite Wand LED will only fit in the Thompson port (extended tip) on the light source turret. Second port is standard ACMI.

NOTE: Standard U.S. power cable included. If a different power cable is required, please contact us at 231-922-0177.

Thompson
SURGICAL INSTRUMENTS



10170 East Cherry Bend Road
Traverse City, Michigan 49684
phone: 231.922.0177
fax: 231.922.0174

thompsonsurgical.com

**Made In
USA**

EC REP

Emergo Europe
Molenstraat 15
2513 BH, The Hague
THE NETHERLANDS

CE 0297

© 2016 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.
Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363,
US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064. Other patents pending.

The Lite Wand LED (40002LWL) and LED Light Source (LLS-2000-2) are manufactured by Sunoptic Technologies® and distributed by Thompson Surgical Instruments.

Rev B
042616
lwlum0416