

Thompson Retractor


Manipulator Handle (#42228)

Versatile and Secure Connection

The Thompson Manipulator Handle attaches to Thompson S-Lock® blades. It provides a versatile and secure connection while allowing the surgeon precise control when positioning the blade. The ergonomic handle provides an ideal angle and grip for the surgeon to comfortably position the blade while maintaining great visibility and exposure.




The Manipulator Handle can be used in addition to a Thompson S-Lock® Retractor Handle to make controlled retraction adjustments to the blade position after the blade has been placed and attached to the frame.



Step 1


With S-Lock retractor handle in locked groove position on blade nipple, connect Manipulator Handle to the top groove on the blade nipple (A).

Once Manipulator Handle is fully over S-Lock nipple, turn blue handle clockwise (B) to tighten to the blade for secure connection.



Step 2

With S-Lock retractor handle detached from frame, or attached to frame in unlocked position (C), retract as desired (D).



Step 3

Once retractor is in place, release Manipulator Handle by turning blue handle counter-clockwise (E) and sliding forward off of S-Lock nipple (F).

 **ORDERING INFO:**
Manipulator Handle #42228


 **NOTICE:**
Manipulator Handle applies only to S-Lock retractor blades.

 **Thompson**
SURGICAL INSTRUMENTS
10341 East Cherry Bend Road
Traverse City, Michigan 49684
phone: 231.922.0177
fax: 231.922.0174
thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2019 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the Thompson "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.
Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064, US9872675. Other patents pending.

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
THE NETHERLANDS

CE0297

Rev C
041819
mhssum

**Made In
USA**