

# Thompson Retractor

## Distractors

### Expose Disc Space

Thompson Surgical now offers Left and Right Caspar Distractors for anterior cervical spine surgery. Left and right versions allow for preferential orientation to be easily attained for left or right-handed surgeons.



*The Caspar Distractor System, the most commonly used distractor system in the world, is considered the gold standard in cervical distraction.*



#### Place

Position distractor in place over screws.



#### Distract

To open distractor, turn the knob.



#### Release

To release tension, press the pawl.



#### NOTE:

The Distractor is a self-retaining device that is used to create space between discs and provide the surgeon a clear view of the vertebral area.



**CERVICAL DISTRACTORS**

REF	ITEM DESCRIPTION	PART #
A	Caspar Distractor Right	743-1000-1
B	Caspar Distractor Left	743-1000-2

**CERVICAL DISTRACTOR PIN KIT**

**#80039**

REF	QTY	ITEM DESCRIPTION	PART #
C	1	Distraction Screwdriver	743-1001-0
D	2	Distraction Screw 12mm	743-0012-0P
	2	Distraction Screw 14mm	743-0014-0P
	2	Distraction Screw 16mm	743-0016-0P

**i NOTE:** Distractors, screws, and screwdriver fit in Instrument Case #50000LW.

**i NOTICE:** Distraction Screws are sold non-sterile. All Distraction Screw sales are final.



Visit us online: [thompsonsurgical.com](http://thompsonsurgical.com)

© 2017 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

® S-Lock®, PLA®, and the "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064

Other patents pending.

Distractors (743-1000-01, 743-1000-02, 743-1011-0), Screws (743-0012-0, 743-0014-0, 743-0016-0), and Screwdriver (743-1001-01) are manufactured by Life Instrument Corporation and distributed by Thompson Surgical Instruments, Inc.

Rev C  
033117  
distsum0417