



Mini
Retractor
Lite LED

Thompson Retractor 



- ADD COOL BRIGHT LED LIGHT TO MOST ANTERIOR CERVICAL BLADES
- ELIMINATE HEAD LAMPS WITH A LOW-PROFILE, REUSABLE LIGHT
- REMOVE SHADOWS FROM OVERHEAD LIGHTS

Uncompromised Exposure.



Mini Retractor Lite LED


Ideal for anterior cervical surgery, the Mini Retractor Lite LED is our lowest profile lighting accessory. Simply attach the included light clip to the light cable and slide on to anterior cervical retractor blades to provide shadowless, cool, LED light in the smallest areas.



MINI RETRACTOR LITE LED KIT (Cervical Use)

QTY	ITEM DESCRIPTION	PART #
	Mini Retractor Lite LED Kit (Cervical Use)	40001MRL
A 1	Mini Retractor Lite LED (Cervical Use)	40002MRL
B 2	Mini Retractor Lite LED - Lite Clip	40026MRL
C 1	Mini Retractor Lite LED Case	50000LED
D 1	LED Light Source	LLS-2000

Exclusive Features



LOW PROFILE LITE CLIP

- Attaches to most anterior cervical blades
- Does not take up room in the incision
- Reusable, 2 included per kit



MINI RETRACTOR LITE LED CABLE

- Low profile design stays out of the way
- Provides cool, bright, LED light
- Reusable



10170 East Cherry Bend Road
 Traverse City, Michigan 49684
 phone: 231.922.0177
 fax: 231.922.0174

thompsonsurgical.com

**For a FREE TRIAL
 call 1.800.227.7543**

Free trial valid for U.S. customers only.
 International customers please call for more information.