


Small Incision Head and Neck System

Lowest Profile for Small Incisions

The Small Incision Head and Neck System is ideal for general and ENT procedures in the head and neck region to expose peripheral vascular structures such as the carotid artery or glands such as the thyroid gland. A wide variety of small retractor blades for skin and shallow tissue retraction are provided for use on a NEW low profile frame designed specifically for the head and neck region. The frame contains two short 11" Elite II Rail Clamps, thinner lateral arms, and our S-Lock Quick Angle Handles with low profile joints.



Frame components in this system have a diameter of $\frac{3}{8}$ " as compared to our standard $\frac{1}{2}$ " frame components.




Low Profile Frame

The Small Incision Head and Neck System includes a NEW low profile and slimmer ($\frac{3}{8}$ " diameter) frame designed specifically for the head and neck region.



Variety of Small Blades

A wide variety of small retractor blades for skin and shallow tissue retraction are provided for Uncompromised Exposure in small incisions.




Quick Angle Handle

Our newest handle, featuring low profile joint and quick angling mechanism. Simply turn the knob with fingers or T-Handle to angle 45° up or down.



SYSTEM NOTICE:

Non-S-Lock system available, but may require a longer lead time. Please call for more information.



S-LOCK SMALL INCISION HEAD AND NECK SYSTEM

#SL90054

REF	QTY	ITEM DESCRIPTION	PART #	REF	QTY	ITEM DESCRIPTION	PART #
A	2	Elite II Rail Clamp with 1 3/8" Cam Joint 10"	41902S	M	1	Volkman Blunt 25mm x 13mm (1" x 1/2")	SL46211
B	2	Angled Arm 19" (12" x 7") 3/8"	44119	N	1	Volkman Sharp 16mm x 13mm (5/8" x 1/2")	SL46212
C	2	Angled Arm 17" (7" x 10") 3/8"	44117	O	1	Volkman Blunt 16mm x 13mm (5/8" x 1/2")	SL46213
D	1	Extension Arm	44014H	P	1	Murphy Sharp 25mm x 13mm (1" x 1/2")	SL46214
E	5	Cam IIP 3/8" Clip-on Quick Angle 8"	SL42127WGP	Q	1	Murphy Blunt 25mm x 13mm (1" x 1/2")	SL46215
F	2	T-Handle	60020	R	1	Murphy Blunt 19mm x 13mm (3/4" x 1/2")	SL46217
G	2	Mini-Richardson 19mm x 25mm (3/4" x 1")	SL46125	S	1	Concave Lateral 25mm x 35mm (1" x 1 1/2")	SL45157A
H	1	Mini-Richardson 19mm x 50mm (3/4" x 2")	SL46125A	T	2	Concave Lateral 25mm x 70mm (1" x 2 3/4")	SL45159A
I	1	Mini-Richardson 25mm x 32mm (1" x 1 1/4")	SL46126	U	2	Concave Medial 25mm x 50mm (1" x 2")	SL45232A
J	2	Richardson 44mm x 51mm (1 3/4" x 2")	SL46149	V	2	Concave Medial 32mm x 70mm (1 1/4" x 2 3/4")	SL45153A
K	1	Richardson 51mm x 102mm (2" x 4")	SL46150		1	Instrument Case 22" x 11" x 3 1/2"	50000G
L	2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115				



Visit us online: thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2017 Thompson Surgical Instruments, Inc. Traverse City, Michigan. Printed in U.S.A.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064
Other patents pending.

Rev A
060917
sp sihns

