

## Liver System Trays

### Keep it Organized

Thompson's Liver Transplant / Resection System contains three instrument trays. This guide outlines which components to place in each tray to remain within the 25lb limit.



### LARGE INSTRUMENT CASE 50000BL

50000BL can contain both Elites and the Rigid Crossbar to equal 22.7 lbs total.

QTY	ITEM DESCRIPTION	PART #
2	Elite II Rail Clamp with 2 Cam Joints 18"	41903AC
1	Bilateral Crossbar	41900F

### STANDARD INSTRUMENT CASE 50000G ( 1 OF 2 )

One 50000G can hold the rest of the frame components including handles, joints, and the angling wrench to equal 22.39 lbs total. ( Pictured on the right. )

QTY	ITEM DESCRIPTION	PART #
1	34" Bilateral Crossbar Hinged (11 1/2" x 8" x 11 1/2")	41900PH
1	18" Angled Arm (8" x 10" @ 45°)	44118N
2	24" Angled Arm	44124W
1	24" Angled Arm (8" x 16")	44124N
2	Cam I Joint Serrated (1/2" x 1/2")	42112C
6	Cam II Clip-on Angling 8"	SL42126WA
2	Micro-Adjustable II Clip-on Angling 10"	SL45006CA
2	Micro-Adjustable II Clip-on Angling 15"	SL45006CAR
2	Angling Head Wrench	43032



### STANDARD INSTRUMENT CASE 50000G ( 2 OF 2 )

The second 50000G has 20lbs (the case weighs 5lbs) to spare to hold the 22 blades in the system (most weigh under a pound).

QTY	ITEM DESCRIPTION	PART #
2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115
2	Malleable 13mm x 152mm (1/2" x 6")	SL46116
2	Malleable 13mm x 203mm (1/2" x 8")	SL46117
1	Malleable 25mm x 203mm (1" x 8")	SL46118
1	Malleable 25mm x 102mm (1" x 4")	SL46118P
1	Malleable 51mm x 203mm (2" x 8")	SL46122
2	Malleable 64mm x 254mm (2 1/2" x 10")	SL46124
1	Malleable 76mm x 254mm (3" x 10")	SO46124A
1	Malleable 102mm x 254mm (4" x 10")	SO46124B
1	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
2	Kelly 89mm x 76mm (3 1/2" x 3")	SO46130
1	Balfour 70mm x 73mm (2 3/4" x 3") [w/small lips]	SO46140C
2	Balfour 83mm x 73mm (3 1/4" x 3")	SO46141
2	Balfour 83mm x 73mm (3 1/4" x 3") [w/small lips]	SO46141C
2	Balfour 83mm x 123mm (3 1/4" x 5")	SO46141D
4	Balfour 32mm x 32mm (1 1/4" x 1 1/4")	SL46142
1	Richardson 57mm x 178mm (2" x 7")	SO46153
1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
1	St. Mark's 64mm x 178mm (2 1/2" x 7")	SO46188
2	Malleable Finger 6"	SL46191
1	Balfour 83mm x 57mm (3 1/4" x 2 1/4")	SO46241A
1	Balfour 100mm x 61mm (4" x 2 3/8")	SO46241B
1	Balfour 114mm x 63mm (4 1/2" x 2 1/2")	SO46241C