


Pediatric Abdominal System

Uncompromised Exposure

Thompson Surgical's Pediatric Abdominal System combines the Infinite Height Rail Clamp, smaller frame arms, and a streamlined selection of small abdominal blades ideal for pediatric patients. Infinite Height Rail Clamp and independent frame arms allow for versatile frame setups and when combined with S-Lock technology blades and handles, the system gives surgeons control over each blade's exact placement while remaining extremely low profile.



*Thompson offers a wide variety of pediatric blade sizes to accommodate your tiniest patients. **

** Custom sizes also available.*




S-Lock® Technology

S-Lock technology reduces unwanted blade swivel and enables controlled, precise blade placement in the smallest, most delicate areas. Easily switch from swivel to locked with the push of a button.



Pediatric Blades

We offer a large selection of blade sizes and styles, catering to a wide range of procedures and patient sizes. Custom sizes also available.




Independent Frame Arms

Frame arms can be independently adjusted to create a low-profile, and provide more control. Independent movement enables safe, versatile retraction.

SYSTEM NOTICE:

Non-S-Lock system still available, but may require a longer lead time. Please call for more information.



PEDIATRIC ABDOMINAL SYSTEM

#SL91034AP

REF	QTY	ITEM DESCRIPTION	PART #
A	1	Infinite Height Rail Clamp with 1 Cam Joint	41902IRC
B	1	Angled Arm 24" (11" x 13")	44124W
C	1	Crossbar with 2 Cam Joints 24"	41900BC
D	2	Angled Arm 12" (6" x 6")	44014S
E	5	Cam II Clip-on Angling 8"	SL42126WA *
F	1	Micro-Adjustable Clip-on Angling 10"	SL45006CA
G	2	Angling Head Wrench	43032
H	2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115
I	2	Malleable 13mm x 152mm (1/2" x 6")	SL46116
J	2	Malleable 13mm x 203mm (1/2" x 8")	SL46117
K	2	Soft Malleable 25mm x 102mm (1" x 4")	SL46118P
L	2	Malleable Pediatric 51mm x 102mm (2" x 4")	SL46122P
M	1	Richardson 19mm x 25mm (3/4" x 1")	SL46125

REF	QTY	ITEM DESCRIPTION	PART #
N	1	Richardson 19mm x 50mm (3/4" x 2")	SL46125A
O	1	Richardson 25mm x 32mm (1" x 1 1/4")	SL46126
P	1	Kelly 38mm x 51mm (1 1/2" x 2")	SL46127
Q	1	Kelly 51mm x 64mm (2" x 2 1/2")	SL46128
R	1	Balfour 76mm x 51mm (3" x 2")	SL46140P
S	2	Balfour 32mm x 32mm (1 1/4" x 1 1/4")	SL46142
T	1	Balfour 51mm x 51mm (2" x 2")	SL46142P
U	1	Richardson 51mm x 102mm (2" x 4")	SL46150
V	1	Harrington 38mm x 127mm (1 1/2" x 5")	SL46162P
W	1	Deaver 25mm x 127mm (1" x 5")	SL46170
	2	Instrument Case w/Pin Mat (22" x 11" x 3 1/2")	50000G
	2	Instrument Case 22" x 11" x 3 1/2" with Pin Mat	50000G

* For infants: Screw Down Joints (#42125N) and Angling Handles (#SL45003AS) recommended.



Visit us online: thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2017 Thompson Surgical Instruments, Inc. Traverse City, Michigan. Printed in U.S.A.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064
Other patents pending.

Rev A
072617
peds

