


# Thompson Retractor

## General / Vascular System

### Uncompromised Exposure

Achieve unsurpassed general and vascular exposure with the Thompson General / Vascular System. The system contains a streamlined selection of abdominal blades ideal for general abdominal and vascular procedures, especially aortic repair. Sterile mount rail clamp and independent frame arms allow for versatile setups and when combined with S-Lock Quick Angle handles, the system gives surgeons control over each blade's exact placement while remaining extremely low profile.



*"Thompson Surgical provides the most hassle-free, versatile retractor. For those complex, difficult procedures, as well as the routine, its ease of setup, strength and use, as well as lack of constant readjustment, always provides me reliable and safe exposure."*

**- Michael Nizzi, DO, FACOS**



#### **S-Lock® Technology**

Reduces unwanted blade swivel, enabling precise blade placement. Switch from swivel to locked with the push of a button.



#### **Quick Angle Technology**

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.




#### **Independent Frame Arms**

Frame arms are independently adjusted, providing more control for safe, versatile, low profile retraction



#### **SYSTEM NOTICE:**

Non-S-Lock system still available, but may require a longer lead time. Please call for more information.



GENERAL / VASCULAR SYSTEM

#SL82003

REF	QTY	ITEM DESCRIPTION	PART #
A	1	Elite III Rail Clamp w/1 Cam Joint 18"	43902AC
B	1	Crossbar w/2 Cam III Joints 24" (11" x 13")	43900BC
C	2	Angled Arm 20" (10" x 10") @ 45°	44120
D	2	Angled Arm 24" (8" x 16") @ 45°	44124N
E	5	Cam II Clip-On Quick Angle 8"	SL42126WGP
F	1	Micro-Adjustable II Clip-On Quick Angle 10"	SL45006CGP
G	2	T-Handle	60020
	1	Instrument Case w/Pin Mat 22" x 11" x 3 1/2"	50000G
	1	Malleable 25mm x 203mm (1" x 8")	SL46118
	2	Conc. Renal Vein 25mm x 140mm (1" x 5 1/2")	SL46119B
H	2	Conc. Renal Vein 25mm x 190mm (1" x 7 1/2")	SL46119C
I	2	Malleable 51mm x 203mm (2" x 8")	SL46122
J	1	Malleable Slotted 51mm x 203mm (2" x 8")	SL46122V
K	1	Malleable 64mm x 254mm (2 1/2" x 10")	SL46124
	1	Malleable 76mm x 254mm (3" x 10")	SO46124A
	1	Malleable Slotted 76mm x 203mm (3" x 8")	SO46124AV
	1	Richardson 19mm x 25mm (3/4" x 1")	SL46125
	2	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
L	1	Kelly 76mm x 89mm (3" x 3 1/2")	SO46130
	1	Fence 102mm x 102mm (4" x 4")	SO46134

REF	QTY	ITEM DESCRIPTION	PART #
	1	Fence 102mm x 127mm (4" x 5")	SO46135
M	1	Fence 102mm x 152mm (4" x 6")	SO46136
	1	Balfour [lips] 70mm x 73mm (2 3/4" x 3")	SO46140
N	1	Balfour [sm lips] 70mm x 73mm (2 3/4" x 3")	SO46140C
	1	Balfour [lips] 83mm x 73mm (3 1/4" x 3")	SO46141
	1	Balfour [sm lips] 83mm x 73mm (3 1/4" x 3")	SO46141C
O	2	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	SL46143
	2	Fenestrated 48mm x 127mm (1 7/8" x 5")	SL46144
P	1	Splanchnich 51mm x 114mm (2" x 4 1/2")	SL46145
	1	Splanchnich 51mm x 152mm (2" x 6")	SL46146
Q	2	Richardson 51mm x 127mm (2" x 5")	SL46151
R	1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
S	1	Deaver 51mm x 127mm (2" x 5")	SL46174
	1	Deaver 76mm x 127mm (3" x 5")	SO46176
T	1	St. Marks [lip] 64mm x 178mm (2 1/2" x 7")	SO46188
U	1	Malleable Finger 152mm (6")	SL46191
	1	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	SO46241A
V	1	Wide Balfour 100mm x 61mm (4" x 2 3/8")	SO46241B
	1	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	SO46241C
	1	Instrument Case w/Pin Mat 22" x 11" x 3 1/2"	50000G



Visit us online: [thompsonsurgical.com](http://thompsonsurgical.com)

For a Free Trial Call **1-800-227-7543** \*

\* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2019 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

® S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc. Patents: tpat.com

Rev C  
101519  
gvss

