


Thompson Retractor

Uncompromised Exposure™



Laparoscopic Selection Guide

Quick and easy set up - Steadier than human hands - Free up staff

Laparoscopic for Elite Rail Clamp

Hands-free

Thompson's laparoscopic products offer hands-free, multi-configuration possibilities for holding scopes and laparoscopic instruments.

"Elite" Laparoscopic products utilize sterile mount Elite or Infinite Height Rail Clamps to mount to OR bed rails.

See pages 6 - 7 for attachment options.



Tip: Use two laparoscopic holding arms with one rail clamp by selecting the Elite III or Infinite Height Rail Clamp with 2 Cam Joints. (shown above with Elite III)



Flexible Holder for Elite or IHRC

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.



Articulating Arm for Elite or IHRC

Holds firmly with a wide range of motion.



Rigid Holder for Elite or IHRC

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.



NOTICE

See pages 6 - 7 for attachments to hold scope or instrument of choice.



ELITE III RAIL CLAMPS

Fast top control attaches to the OR table in the sterile field. Additional heights available, call for information.

REF	ITEM DESCRIPTION	PART #
A	18" (45cm) Tall with 1 Cam Joint	43902AC
	18" (45cm) Tall with 2 Cam Joints	43903AC

INFINITE HEIGHT RAIL CLAMPS (IHRC)

The Infinite Height Rail Clamp adjusts between 13" - 22" tall to accommodate all patient sizes and remain low profile.

REF	ITEM DESCRIPTION	PART #
B	IHRC with 1 Cam Joint	43902IRC
	IHRC with 2 Cam Joints	43903IRC



LAPAROSCOPIC ARTICULATING ARM HOLDER

Holds firmly with a wide range of motion.

ITEM DESCRIPTION	PART #
Laparoscopic Articulating Arm for Elite (or IHRC)	42138LE



FLEXIBLE LAPAROSCOPIC HOLDER

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.

ITEM DESCRIPTION	PART #
Flexible Laparoscopic Holder for Elite (or IHRC)	42133DEN



RIGID LAPAROSCOPIC HOLDER

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.

ITEM DESCRIPTION	PART #
Universal Rigid Laparoscopic Holder for Elite (or IHRC)	42139EU

NOTICE
Elite II Rail Clamps are still available, call for information.

NOTICE
Rigid Laparoscopic Holder with integrated Duo Head or S-Lock® Nathanson Adapter also available, call for information.

Laparoscopic for Power Rail Clamp

Lower profile, basic design

Thompson's original laparoscopic holding system, the Power Rail Clamp set up, is still our lowest profile laparoscopic holding option. It holds securely and hands-free, while remaining low profile and out of the way.

See pages 6 - 7 for attachment options.



Tip: Use two Power Rail Clamps with laparoscopic holding arms to hold scopes and instruments simultaneously. (shown above)



Flexible Holder for PRC

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.



Articulating Arm for PRC

Holds firmly with a wide range of motion.



Rigid Holder for PRC

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.



NOTICE

See pages 6 - 7 for attachments to hold scope or instrument of choice.



POWER RAIL CLAMP

The Power Rail Clamp is ideal for those who prefer a lower profile, basic rail attachment.

ITEM DESCRIPTION	PART #
Power Rail Clamp (PRC)	41950



FLEXIBLE LAPAROSCOPIC HOLDER

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.

ITEM DESCRIPTION	PART #
Flexible Laparoscopic Holder for Power Rail Clamp	42133DN



LAPAROSCOPIC ARTICULATING ARM HOLDER

Holds firmly with a wide range of motion.

ITEM DESCRIPTION	PART #
Laparoscopic Articulating Arm for Power Rail Clamp	42138L



RIGID LAPAROSCOPIC HOLDER

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.

ITEM DESCRIPTION	PART #
Universal Rigid Laparoscopic Holder for PRC	42139U



NOTICE

Rigid Laparoscopic Holder with integrated Duo Head or S-Lock® Nathanson Adapter also available, call for information.

Laparoscopic Attachments

Attach your scope or instrument

A variety of attachments connect quickly and easily, while protecting scopes.

Thompson's laparoscopic attachments easily and securely connect to any of the holding arms shown on pages 2 - 5.




"Hands down, this should be every surgeon's go-to retractor."

- Michael Nizzi, DO, FACOS



Adjustable Scope Grip

- Easily adjusts to fit 5-11mm scopes
- Holds scope firmly in position, while protecting scope shaft
- When loosened, allows for panning



Scope Clip

- Allows quick in / out panning of flat scopes
- Clips may be custom ordered to any size



S-Lock® Nathanson Hooks

- Eliminates slippage and rotation
- Holds securely
- Quick and easy adjustments



NOTICE

Duo Head and Scope Grip attachments also available, see page 7 for information.



ADJUSTABLE SCOPE GRIP

Holds scope firmly while protecting shaft. Fits scopes 5 - 11mm.

ITEM DESCRIPTION	PART #
Adjustable Scope Grip 5 - 11mm	33181



SCOPE GRIP AND BUSHINGS

Holds scope tightly while preventing rotation of angled scopes.

REF	ITEM DESCRIPTION	PART #
A	Scope Grip	33154
B	Scope Grip Bushing 5mm **	42136F
	Scope Grip Bushing 10mm	42136T


** Additional sizes from 5.1mm - 5.8mm also available



DUO HEAD

Holds 3mm, 5mm, or 10mm instruments firmly in place.

ITEM DESCRIPTION	PART #
Duo Head with Tool Grasp Knob 3/5/10mm	42160N



SCOPE CLIP *

Quick in and out panning of flat scopes. Custom sizes available.

REF	ITEM DESCRIPTION	PART #
C	Scope Clip 5mm	42137F
	Scope Clip 7mm	42137S
	Scope Clip 10mm	42137T
	Scope Clip 11mm	42137E

* Scope Clip may be custom ordered to fit any size scope diameter



S-LOCK® NATHANSON HOOKS

Eliminates slippage and rotation, quick and easy adjustment.

REF	ITEM DESCRIPTION	PART #
	Small 5mm (169mm x 125mm x 12mm)	SL42200
D	Medium 5mm (203mm x 167mm x 16mm)	SL42201
	Large 5mm (240mm x 222mm x 21mm)	SL42202
	Small 6 ½mm (169mm x 125mm x 12mm)	SL42203
	Medium 6 ½mm (203mm x 167mm x 16mm)	SL42204
	Large 6 ½mm (240mm x 222mm x 21mm)	SL42205
E	S-Lock Nathanson Adapter (required for S-Lock Nathanson Hook use)	SL42161

For a Free Trial Call Today*

1.800.227.7543

Rev G
021418
lapsg



Uncompromised Exposure

VISIT US ONLINE: THOMPSONSURGICAL.COM

Thompson
SURGICAL INSTRUMENTS

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064
Other patents pending.

Made In
USA