

# Thompson Retractor

Uncompromised Exposure™



## Laparoscopic Selection Guide

Quick and easy set up - Steadier than human hands - Free up staff

## Laparoscopic for Elite Rail Clamp

### Hands-free

Thompson's laparoscopic products offer hands-free, multi-configuration possibilities for holding scopes and laparoscopic instruments.

"Elite" Laparoscopic products utilize sterile mount Elite or Infinite Height Rail Clamps to mount to OR bed rails.

See pages 6 - 7 for attachment options.



**Tip:** Use two laparoscopic holding arms with one rail clamp by selecting the Elite III or Infinite Height Rail Clamp with 2 Cam Joints. ( shown above with Elite III )



#### Flexible Holder for Elite or IHRC

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.



#### Articulating Arm for Elite or IHRC

Holds firmly with a wide range of motion.



#### Rigid Holder for Elite or IHRC

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.



#### NOTICE

See pages 6 - 7 for attachments to hold scope or instrument of choice.



### ELITE III RAIL CLAMPS

Fast top control attaches to the OR table in the sterile field. Additional heights available, call for information.

REF	ITEM DESCRIPTION	PART #
A	18" (45cm) Tall with 1 Cam Joint	43902AC
	18" (45cm) Tall with 2 Cam Joints	43903AC

### INFINITE HEIGHT RAIL CLAMPS ( IHRC )

The Infinite Height Rail Clamp adjusts between 13" - 22" tall to accommodate all patient sizes and remain low profile.

REF	ITEM DESCRIPTION	PART #
B	IHRC with 1 Cam Joint	43902IRC
	IHRC with 2 Cam Joints	43903IRC



### LAPAROSCOPIC ARTICULATING ARM HOLDER

Holds firmly with a wide range of motion.

ITEM DESCRIPTION	PART #
Laparoscopic Articulating Arm for Elite ( or IHRC )	42138LE



### FLEXIBLE LAPAROSCOPIC HOLDER

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.

ITEM DESCRIPTION	PART #
Flexible Laparoscopic Holder for Elite ( or IHRC )	42133DEN



### RIGID LAPAROSCOPIC HOLDER

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.

ITEM DESCRIPTION	PART #
Universal Rigid Laparoscopic Holder for Elite ( or IHRC )	42139EU

**NOTICE**  
Elite II Rail Clamps are still available, call for information.

**NOTICE**  
Rigid Laparoscopic Holder with integrated Duo Head or S-Lock® Nathanson Adapter also available, call for information.

## Laparoscopic for Power Rail Clamp

### Lower profile, basic design

Thompson's original laparoscopic holding system, the Power Rail Clamp set up, is still our lowest profile laparoscopic holding option. It holds securely and hands-free, while remaining low profile and out of the way.

See pages 6 - 7 for attachment options.



*Tip: Use two Power Rail Clamps with laparoscopic holding arms to hold scopes and instruments simultaneously. ( shown above )*



#### Flexible Holder for PRC

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.



#### Articulating Arm for PRC

Holds firmly with a wide range of motion.



#### Rigid Holder for PRC

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.



#### NOTICE

See pages 6 - 7 for attachments to hold scope or instrument of choice.



### POWER RAIL CLAMP

The Power Rail Clamp is ideal for those who prefer a lower profile, basic rail attachment.

ITEM DESCRIPTION	PART #
Power Rail Clamp ( PRC )	41950



### FLEXIBLE LAPAROSCOPIC HOLDER

Unlimited range of motion and positioning with adjustable tension from flexible to rigid.

ITEM DESCRIPTION	PART #
Flexible Laparoscopic Holder for Power Rail Clamp	42133DN



### LAPAROSCOPIC ARTICULATING ARM HOLDER

Holds firmly with a wide range of motion.

ITEM DESCRIPTION	PART #
Laparoscopic Articulating Arm for Power Rail Clamp	42138L



### RIGID LAPAROSCOPIC HOLDER

Precise, rigid positioning for instruments requiring few adjustments or holding excess weight.

ITEM DESCRIPTION	PART #
Universal Rigid Laparoscopic Holder for PRC	42139U



#### NOTICE

Rigid Laparoscopic Holder with integrated Duo Head or S-Lock® Nathanson Adapter also available, call for information.

## Laparoscopic Attachments

### Attach your scope or instrument

A variety of attachments connect quickly and easily, while protecting scopes.

Thompson's laparoscopic attachments easily and securely connect to any of the holding arms shown on pages 2 - 5.



*"Hands down, this should be every surgeon's go-to retractor."*

- Michael Nizzi, DO, FACOS



#### Adjustable Scope Grip

- Easily adjusts to fit 5-6mm or 10-11mm scopes
- Holds scope firmly in position, while protecting scope shaft
- When loosened, allows for panning



#### Scope Clip

- Allows quick in / out panning of 0° scopes
- Scope Clip may be ordered to fit any size scope diameter



#### S-Lock® Nathanson Hooks

- Eliminates slippage and rotation
- Holds securely
- Quick and easy adjustments



#### NOTICE

Duo Head and Scope Grip attachments also available, see page 7 for information.





### ADJUSTABLE SCOPE GRIP

Holds scope firmly while protecting shaft. Fits scopes 5-6mm or 10-11mm.

ITEM DESCRIPTION	PART #
Adjustable Scope Grip 5 - 11mm	33181



### SCOPE GRIP AND BUSHINGS

Holds scope tightly while preventing rotation of angled scopes.

REF	ITEM DESCRIPTION	PART #
A	Scope Grip	33154
B	Scope Grip Bushing 5mm **	42136F
	Scope Grip Bushing 10mm	42136T


\*\* Additional sizes from 5.1mm - 5.8mm also available



### DUO HEAD

Holds 3mm, 5mm, or 10mm instruments firmly in place.

ITEM DESCRIPTION	PART #
Duo Head with Tool Grasp Knob 3/5/10mm	42160N



### SCOPE CLIP \*

Quick in and out panning of 0° scopes. Reusable.

REF	ITEM DESCRIPTION	PART #
C	Scope Clip 5mm	42137F
	Scope Clip 7mm	42137S
	Scope Clip 10mm	42137T
	Scope Clip 11mm	42137E

\*Scope Clip may be ordered to fit any size scope diameter



### S-LOCK® NATHANSON HOOKS

Eliminates slippage and rotation, quick and easy adjustment.

REF	ITEM DESCRIPTION	PART #
	Small 5mm (169mm x 125mm x 12mm)	SL42200
D	Medium 5mm (203mm x 167mm x 16mm)	SL42201
	Large 5mm (240mm x 222mm x 21mm)	SL42202
	Small 6 ½mm (169mm x 125mm x 12mm)	SL42203
	Medium 6 ½mm (203mm x 167mm x 16mm)	SL42204
	Large 6 ½mm (240mm x 222mm x 21mm)	SL42205
E	S-Lock Nathanson Adapter ( required for S-Lock Nathanson Hook use )	SL42161

For a Free Trial Call Today\*

**1.800.227.7543**

Rev H  
121021  
lapsg




Uncompromised Exposure

VISIT US ONLINE: [THOMPSONSURGICAL.COM](http://THOMPSONSURGICAL.COM)

**Thompson**   
SURGICAL INSTRUMENTS

\* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2021 Thompson Surgical Instruments, Inc.

® S-Lock®, PLA®, and the "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.  
Patents: [tapat.com](http://tapat.com)

**Made In**  
**USA** 