


# Thompson Retractor

## Pediatric System

### Uncompromised Exposure

The Thompson Pediatric System combines the Infinite Height Rail Clamp, smaller frame arms, and a streamlined selection of small abdominal blades ideal for pediatric patients. The combination of versatile frame setups and S-Lock technology blades and handles gives surgeons control over each blade's exact placement while remaining extremely low profile.



*"The Thompson Retractor makes liver transplant easy. It provides outstanding exposure for pediatric liver transplant, the best retractor I have seen in my 25 years as a transplant surgeon."*

**- Michael Wachs, MD**



### S-Lock® Technology

Reduces unwanted blade swivel, enabling precise blade placement. Switch from swivel to locked with the push of a button.



### Pediatric Blades

We offer a wide variety of pediatric blade sizes to accommodate your tiniest patients. Custom sizes also available.




### Independent Frame Arms

Frame arms can be independently adjusted to create a low-profile, and provide more control. Independent movement enables safe, versatile retraction.



#### SYSTEM NOTICE:

Non-S-Lock system still available, but may require a longer lead time. Please call for more information.



PEDIATRIC SYSTEM

#SL82006

REF	QTY	ITEM DESCRIPTION	PART #
A	1	Infinite Height Rail Clamp w/1 Cam III Joint	43902IRC
B	1	Crossbar with 2 Cam III Joints 24"	43900BC
C	2	Angled Arm 20" (10" x 10")	44120
D	2	Angled Arm 12" (6" x 6")	44014S
E	5	Cam II Clip-on Quick Angle 8" *	SL42126WGP
F	1	Micro-Adjustable Clip-on Quick Angle 10"	SL45006CGP
G	2	T-Handle	60020
H	2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115
I	2	Malleable 13mm x 152mm (1/2" x 6")	SL46116
J	2	Malleable 13mm x 203mm (1/2" x 8")	SL46117
K	2	Soft Malleable 25mm x 102mm (1" x 4")	SL46118P
L	2	Malleable 51mm x 102mm (2" x 4")	SL46122P
M	1	Malleable Finger 76mm (3")	SL46191SP

REF	QTY	ITEM DESCRIPTION	PART #
N	1	Malleable Finger 102mm (4")	SL46191P
O	1	Richardson 19mm x 25mm (3/4" x 1")	SL46125
P	1	Richardson 19mm x 50mm (3/4" x 2")	SL46125A
Q	1	Richardson 25mm x 32mm (1" x 1 1/4")	SL46126
R	1	Kelly 38mm x 51mm (1 1/2" x 2")	SL46127
S	1	Kelly 51mm x 64mm (2" x 2 1/2")	SL46128
T	1	Balfour 76mm x 51mm (3" x 2")	SL46140P
U	2	Balfour 32mm x 32mm (1 1/4" x 1 1/4")	SL46142
V	1	Balfour 51mm x 54mm (2" x 2")	SL46142P
W	1	Richardson 51mm x 102mm (2" x 4")	SL46150
X	1	Harrington 38mm x 127mm (1 1/2" x 5")	SL46162P
Y	1	Deaver 25mm x 127mm (1" x 5")	SL46170
	2	Instrument Case w/Pin Mat (22" x 11" x 3 1/2")	50000G

\* For infants: Screw Down Joints (#42125N) and Quick Angle Handles (#SL45003GS) recommended.



Visit us online: [thompsonsurgical.com](http://thompsonsurgical.com)

For a Free Trial Call **1-800-227-7543** \*

\* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064  
Other patents pending.

Rev A  
031318  
pedss

