

Small Incision System

Lowest Profile for Head and Neck Procedures

The Small Incision System is ideal for general and ENT procedures in the head and neck region to expose peripheral vascular structures including the carotid artery, or glands such as the thyroid gland.




"Thompson products always facilitated my operations!! Great engineering and quality!"

- John (Jeb) Hallett, MD



Variety of Small Blades

System includes a large selection of small retractor blades for skin and shallow tissue retraction.



NEW Low Profile Frame

Slim, low profile frame is designed specifically for the head and neck region, and includes slimmer ($\frac{3}{8}$ " diameter) arms and low profile handles.



Quick Angle Handles

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.



SYSTEM NOTICE:

Frame components in this system have a diameter of $\frac{3}{8}$ " as compared to our standard $\frac{1}{2}$ " frame components.



SMALL INCISION SYSTEM

#SL82009

REF	QTY	ITEM DESCRIPTION	PART #
A	2	Elite III Rail Clamp w/1 (3/8") Cam Joint 11"	43902S
B	2	19" Angled Arm (12" x 7") @ 45° (3/8")	44119
C	2	17" Angled Arm (10" x 7") @ 45° (3/8")	44117
D	1	3/8" Extension Arm 12" with Hinge Joint	44014H
E	5	Cam II Clip-on Quick Angle 8" (3/8")	SL42127WGP
F	2	T-Handle	60020
G	2	Richardson 19mm x 25mm (3/4" x 1")	SL46125
H	1	Richardson 19mm x 51mm (3/4" x 2")	SL46125A
I	1	Richardson 25mm x 32mm (1" x 1 1/4")	SL46126
J	2	Richardson 44mm x 51mm (1 3/4" x 2")	SL46149
K	1	Richardson 51mm x 102mm (2" x 4")	SL46150
L	2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115

REF	QTY	ITEM DESCRIPTION	PART #
M	1	Volkman Blunt 25mm x 13mm (1" x 1/2")	SL46211
N	1	Volkman Sharp 16mm x 13mm (5/8" x 1/2")	SL46212
O	1	Volkman Blunt 16mm x 13 mm (5/8" x 1/2")	SL46213
P	1	Murphy Sharp 25mm x 13mm (1" x 1/2")	SL46214
Q	1	Murphy Blunt 25mm x 13mm (1" x 1/2")	SL46215
R	1	Murphy Blunt 19mm x 13mm (3/4" x 1/2")	SL46217
S	1	Concave Lateral 25mm x 40mm (1" x 1 1/2")	SL45157A
T	2	Concave Lateral 25mm x 70mm (1" x 2 3/4")	SL45159A
U	2	Concave Medial 25mm x 50mm (1" x 2")	SL45232A
V	2	Concave Medial 32mm x 70mm (1 1/4" x 2 3/4")	SL45153A
	1	Instrument Case	50000G



Visit us online: thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064

Other patents pending.

Rev A
010818
siss

