

## Balfour Blade Styles

### Find the Style for Your Procedure

Thompson Surgical offers many sizes and styles of Balfour blades. Refer to the images below in order to select the style appropriate to your application.



### BALFOUR BLADES - STANDARD LIPS


ITEM DESCRIPTION	PART #
32mm x 32mm (1 1/4" x 1 1/4")	SL46142
51mm x 54mm (2" x 2")	SL46142P
76mm x 51mm (3" x 2")	SL46140P
70mm x 73mm (2 3/4" x 3")	SO46140
83mm x 73mm (3 1/4" x 3")	SO46141
83mm x 123mm (3 1/4" x 5")	SO46141D
83mm x 175mm (3 1/4" x 6 7/8")	SO46141F
<b>Radiolucent</b> 65mm x 72mm (2 1/2" x 2 3/4")	SL46140ET *
<b>Radiolucent</b> 65mm x 72mm (2 1/2" x 2 3/4")	SO46140ET
<b>Radiolucent</b> 83mm x 70mm (3 1/4" x 2 3/4")	SL46141ET *
<b>Radiolucent</b> 83mm x 70mm (3 1/4" x 2 3/4")	SO46141ET
<b>Radiolucent</b> 65mm x 72mm (2 1/2" x 2 3/4") Anodized Aluminum	SO46140EB
<b>Radiolucent</b> 83mm x 70mm (3 1/4" x 2 3/4") Anodized Aluminum	SO46141EB
<b>Radiolucent</b> 32mm x 32mm (1 1/4" x 1 1/4") Anodized Aluminum	SL46142EB
<b>Radiolucent</b> 83mm x 127mm (3 1/4" x 5") Anodized Aluminum	SO46141DEB

\* For anterior lumbar use only.

**i NOTICE:** Any Balfour blades may be ordered with small lips or without lips. Call to order.

### BALFOUR BLADES - SMALL LIPS

ITEM DESCRIPTION	PART #
70mm x 73mm (2 3/4" x 3") [with small lips]	SO46140C
83mm x 73mm (3 1/4" x 3") [with small lips]	SO46141C
<b>Radiolucent</b> 65mm x 72mm (2 1/2" x 2 3/4") [with small lips]	SO46140CET
<b>Radiolucent</b> 83mm x 70mm (3 1/4" x 2 3/4") [with small lips]	SO46141CET
<b>Radiolucent</b> 65mm x 72mm (2 1/2" x 2 3/4") [w/small lips] Anodized Aluminum	SO46140CEB
<b>Radiolucent</b> 83mm x 70mm (3 1/4" x 2 3/4") [w/small lips] Anodized Aluminum	SO46141CEB



### WIDE BALFOUR ( ROCHARD ) BLADES - NO LIPS

ITEM DESCRIPTION	PART #
83mm x 57mm (3 1/4" x 2 1/4")	SO46241A
100mm x 61mm (4" x 2 3/8")	SO46241B
114mm x 63mm (4 1/2" x 2 1/2")	SO46241C
<b>Radiolucent</b> 83mm x 57mm (3 1/4" x 2 1/4") Anodized Aluminum	SO46241AEB
<b>Radiolucent</b> 102mm x 57mm (4" x 2 1/4") Anodized Aluminum	SO46241BEB
<b>Radiolucent</b> 114mm x 57mm (4 1/2" x 2 1/4") Anodized Aluminum	SO46241CEB