


Thompson Retractor

GYN Oncology System

Uncompromised Exposure for GYN Oncology

The Thompson GYN Oncology System provides Uncompromised Exposure for the most complex and widespread gynecological oncology procedures. System can also be adapted for simple abdominal and vaginal procedures.




Versatile system for use in gynecological oncology or simple GYN procedures.



Quick Angle Handles

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.



Bilateral Stability and Rigidity

Stability for wide exposures and asymmetrical costal margin retraction.




S-Lock® Technology

S-Lock technology reduces unwanted blade swivel and enables precise blade placement. Switch from swivel to locked with the push of a button.



SYSTEM NOTICE:

Smaller, unilateral GYN System (#SL82010) also available.



GYN ONCOLOGY SYSTEM

#SL82011

REF	QTY	ITEM DESCRIPTION	PART #
A	2	Elite III Rail Clamp with 2 Cam Joints 18"	43903AC
B	1	34" Bilateral Crossbar Hinged (11 1/2" x 8 1/2" x 11 1/2")	41900PH
C	1	Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4"	41900F
D	1	Crossbar w/2 Cam Joints 24"	43900BC
E	2	24" Angled Arm (11" x 13") @ 45°	44124W
F	2	18" Angled Arm (8" x 10") @ 45°	44118N
G	1	Lower Abdominal Bar 14" x 6"	44112SP
H	6	Cam II Clip-on Quick Angle 8"	SL42126WGP
I	2	Micro-Adjustable II Clip-on Quick Angle 10"	SL45006CGP
	2	Micro-Adjustable II Clip-on Quick Angle 15"	SL45006CGRP
J	2	T-Handle	60020
	1	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	SO46241A
K	1	Wide Balfour 100mm x 61mm (4" x 2 3/8")	SO46241B
	1	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	SO46241C
L	1	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	SO46140C
	2	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	SO46141D
	2	Balfour [small lips] 83mm x 73mm (3 1/4" x 3")	SO46141C
M	1	Deaver 51mm x 127mm (2" x 5")	SL46174

REF	QTY	ITEM DESCRIPTION	PART #
N	1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
	1	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
O	2	Kelly 89mm x 76mm (3 1/2" x 3")	SO46130
	1	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115
P	1	Malleable 13mm x 152mm (1/2" x 6")	SL46116
	2	Malleable 13mm x 203mm (1/2" x 8")	SL46117
	1	Soft Malleable 25mm x 102mm (1" x 4")	SL46118P
Q	1	Malleable 25mm x 203mm (1" x 8")	SL46118
	1	Malleable 51mm x 203mm (2" x 8")	SL46122
	2	Malleable 64mm x 254mm (2 1/2" x 10")	SL46124
	1	Malleable 76mm x 254mm (3" x 10")	SO46124A
R	1	Malleable 102mm x 254mm (4" x 10")	SO46124B
S	2	Malleable Finger 152mm (6")	SL46191
T	2	Richardson 51mm x 127mm (2" x 5")	SO46151
	1	Richardson 57mm x 178mm (2" x 7")	SO46153
U	1	St. Mark's 64mm x 178mm (2 1/2" x 7")	SO46188
	2	Instrument Case w/Pin Mat (22" x 11" x 3 1/2")	50000G
	1	Instrument Case (26" x 10" x 5")	50000BL



Visit us online: thompsonsurgical.com · For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064
Other patents pending.

Rev A
010918
gynoncss

