


Thompson Retractor

GYN Oncology System

Uncompromised Exposure for GYN Oncology

The Thompson GYN Oncology System provides Uncompromised Exposure for the most complex and widespread gynecological oncology procedures. System can also be adapted for simple abdominal and vaginal procedures.




Versatile system for use in gynecological oncology or simple GYN procedures.



Quick Angle Handles

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.



Bilateral Stability and Rigidity

Stability for wide exposures and asymmetrical costal margin retraction.



S-Lock® Technology

S-Lock technology reduces unwanted blade swivel and enables precise blade placement. Switch from swivel to locked with the push of a button.



SYSTEM NOTICE:

Smaller, unilateral GYN System (#SL82010) also available.

GYN ONCOLOGY SYSTEM

#SL82011

| QTY | ITEM DESCRIPTION | PART # |
|-----|---|-------------|
| 2 | Elite III Rail Clamp with 2 Cam Joints 18" | 43903AC |
| 1 | Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4" | 41900F |
| 1 | Lower Abdominal Bar 20" (14" x 6") | 44112SP |
| 1 | Crossbar w/2 Cam Joints 24" | 43900BC |
| 1 | Instrument Case (26" x 10" x 5") | 50000BL |
| 1 | Bilateral Crossbar Hinged 34" (11 1/2" x 8 1/2" x 11 1/2") | 41900PH |
| 1 | Cam III Serrated Slide-On Joint | 43112C |
| 2 | Angled Arm 18" (8" x 10") @ 45° | 44118N |
| 2 | Angled Arm 24" (11" x 13") @ 45° | 44124W |
| 6 | Cam II Clip-on Quick Angle 8" | SL42126WGP |
| 2 | Micro-Adjustable II Clip-on Quick Angle 10" | SL45006CGP |
| 2 | Micro-Adjustable II Clip-on Quick Angle 15" | SL45006CGRP |
| 2 | T-Handle | 60020 |
| 1 | Instrument Case (23" x 11" x 3 1/2") | 50000G |
| 1 | Malleable 13mm x 64mm (1/2" x 2 1/2") | SL46115 |
| 1 | Malleable 13mm x 152mm (1/2" x 6") | SL46116 |
| 2 | Malleable 13mm x 203mm (1/2" x 8") | SL46117 |
| 1 | Malleable 25mm x 203mm (1" x 8") | SL46118 |
| 1 | Malleable 25mm x 254mm (1" x 10") | SL46118A |
| 1 | Malleable 25mm x 102mm (1" x 4") | SL46118P |
| 1 | Malleable 51mm x 203mm (2" x 8") | SL46122 |
| 2 | Malleable 64mm x 254mm (2 1/2" x 10") | SL46124 |
| 1 | Malleable 76mm x 254mm (3" x 10") | SO46124A |
| 1 | Malleable 102mm x 254mm (4" x 10") | SO46124B |
| 1 | Kelly 38mm x 51mm (1 1/2" x 2") | SL46127 |


| QTY | ITEM DESCRIPTION | PART # |
|-----|---|----------|
| 1 | Kelly 51mm x 64mm (2" x 2 1/2") | SL46128 |
| 1 | Kelly 64mm x 76mm (2 1/2" x 3") | SO46129 |
| 2 | Kelly 76mm x 89mm (3 1/2" x 3") | SO46130 |
| 1 | Fence 102mm x 127mm (4" x 5") | SO46135 |
| 1 | Balfour [small lips] 70mm x 73mm (2 3/4" x 3") | SO46140C |
| 2 | Balfour [small lips] 83mm x 73mm (3 1/4" x 3") | SO46141C |
| 2 | Balfour [small lips] 83mm x 123mm (3 1/4" x 5") | SO46141D |
| 2 | Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4") | SL46142 |
| 2 | Fenestrated 44mm x 89mm (1 3/4" x 3 1/2") | SL46143 |
| 1 | Richardson 51mm x 102mm (2" x 4") | SL46150 |
| 1 | Richardson 51mm x 127mm (2" x 5") | SO46151 |
| 1 | Richardson 57mm x 178mm (2" x 7") | SO46153 |
| 1 | Harrington 64mm x 152mm (2 1/2" x 6") | SL46160 |
| 1 | Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2") | SO46168B |
| 1 | Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2") | SO46169B |
| 1 | Deaver 51mm x 127mm (2" x 5") | SL46174 |
| 1 | Heany 25mm x 102mm (1" x 4") | SL46180 |
| 1 | St. Mark's 64mm x 178mm (2 1/2" x 7") | SO46188 |
| 2 | Malleable Finger 152mm (6") | SL46191 |
| 1 | Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4") | SO46241A |
| 1 | Wide Balfour 100mm x 61mm (4" x 2 3/8") | SO46241B |
| 1 | Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2") | SO46241C |
| 2 | Double Angled 25mm x 76mm (1" x 3") | SL46660 |
| 2 | Double Angled 13mm x 70mm (1/2" x 2 3/4") | SL46663 |
| 1 | Instrument Case (23" x 11" x 3 1/2") | 50000G |



Visit us online: thompsonsurgical.com · For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2020 Thompson Surgical Instruments, Inc.

® S-Lock®, PLA®, and the "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc. Patents: tpat.com

