

Liver / Oncology System Trays

Keep it Organized

Thompson's Liver / Oncology System (#SL82007) * contains three instrument trays. This guide outlines which components to place in each tray to remain within the 25 lb limit.



LARGE INSTRUMENT CASE 50000BL

50000BL can contain both Elites, bilateral crossbar, and Lower Abdominal Bar to equal 22.8 lbs total.

QTY	ITEM DESCRIPTION	PART #
2	Elite III Rail Clamp with 2 Joints 18"	43903AC
1	34" Bilateral Crossbar Hinged (11 1/2" x 8" x 11 1/2")	41900PH
1	Lower Abdominal Bar 20" (14" x 6")	44112SP

STANDARD INSTRUMENT CASE 50000G (1 OF 2)

One 50000G can hold the rest of the frame components including handles, joints, and T-Handles to equal 23 lbs total. (Pictured above right.)

QTY	ITEM DESCRIPTION	PART #
1	18" Angled Arm (8" x 10" @ 45°)	44118N
2	24" Angled Arm (11" x 13" @ 45°)	44124W
2	Cam Joint Serrated (1/2" x 1/2")	43112C
6	Cam II Clip-on Quick Angle 8"	SL42126WGP
2	Micro-Adjustable II Clip-on Quick Angle 10"	SL45006CGP
2	Micro-Adjustable II Clip-on Quick Angle 15"	SL45006CGRP
2	T-Handle	60020



STANDARD INSTRUMENT CASE 50000G (2 OF 2)

The second 50000G holds the 26 blades in the system, totaling 14.4 lbs.

QTY	ITEM DESCRIPTION	PART #
1	Malleable 51mm x 203mm (2" x 8")	SL46122
2	Malleable 64mm x 254mm (2 1/2" x 10")	SL46124
1	Malleable 76mm x 254mm (3" x 10")	SO46124A
1	Malleable 102mm x 254mm (4" x 10")	SO46124B
1	Kelly 38mm x 51mm (1 1/2" x 2")	SL46127
1	Kelly 51mm x 64mm (2" x 2 1/2")	SL46128
1	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
1	Kelly 76mm x 89mm (3" x 3 1/2")	SO46130
2	Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	SL46142
1	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	SO46140C
2	Balfour [small lips] 83mm x 73mm (3 1/4" x 3")	SO46141C
2	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	SO46141D
1	Richardson 51mm x 127mm (2" x 5")	SO46151
1	Richardson 51mm x 178mm (2" x 7")	SO46153
1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
1	St Marks [with lip] 64mm x 178mm (2 1/2" x 7")	SO46188
1	Malleable Finger 152mm (6")	SL46191
1	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	SO46241A
1	Wide Balfour 100mm x 61mm (4" x 2 3/8")	SO46241B
1	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	SO46241C
1	Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2")	SO46169B
1	Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2")	SO46168B

* Guidelines also apply to Radiolucent Liver / Oncology System (#SL82008).