


# Thompson Retractor

## Handles for Non-Thompson Systems

### Adapt any Thompson® blade to your frame

Thompson offers a wide range of innovative features providing Uncompromised Exposure™.

Without upgrading the entire system, users of other retractors may incorporate some of this technology with the use of Thompson Adapter Handles.



### Receive a Discount When You Trade In Your Retractor

*Try any of our new components with a FREE TRIAL. If you love them, you can trade in your existing components for a generous discount.*



#### S-Lock® Technology

Reduces unwanted blade swivel, enabling precise blade placement. Switch from swivel to locked with the push of a button.




#### Quick Angle Technology

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.



#### Wide Variety of Blades


Thompson offers a large selection of both radiolucent and radiopaque blades, catering to multiple procedures and patient sizes. Custom blades also available, call for information.



**RATCHET RETRACTOR HANDLES \***


REF	ITEM DESCRIPTION	PART #
A	Quick Angle 10 1/4"	SL45005GL
B	10 1/4" (non-angling)	SL45005L
C	Angling 10 1/4"	SL45005AL

\* Compatible with Bookwalter® brand retractors.



**THOMPSON ADAPTER HANDLES FOR SYNFRAME®**

REF	ITEM DESCRIPTION	PART #
D	Quick Angle 8 3/4"	SL45014GL
E	8 3/4" (non-angling)	SL45014L
F	Angling 8 3/4"	SL45014AL



**THOMPSON ADAPTER HANDLES FOR OMNI-TRACT®**

REF	ITEM DESCRIPTION	PART #
G	Quick Angle 8 3/4"	SL45004GL
H	8 3/4" (non-angling)	SL45004L
I	S-Lock Angling 8 3/4"	SL45004AL

**NOTICE:**

Bookwalter® is a Registered Trademark of Symmetry Surgical®. Thompson Surgical Instruments, Inc. has no affiliation with the Bookwalter® Retractor or Symmetry Surgical®.

SynFrame® is a Registered Trademark of DePuy Synthes, Inc®. Thompson Surgical Instruments, Inc. has no affiliation with the SynFrame® or DePuy Synthes, Inc®.

Omni-Tract® is a Registered Trademark of Integra®. Thompson Surgical Instruments, Inc. has no affiliation with Omni-Tract® or Integra®.



Visit us online: [thompsonsurgical.com](http://thompsonsurgical.com)

For a Free Trial Call **1-800-227-7543** \*

\* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan.

© S-Lock®, PLA®, Thompson® and the "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064  
Other patents pending.

Rev A  
040518  
compss

