


Thompson Retractor

Trauma System


Fast, Easy Set Up

Sterile mount rail clamp and independent frame arms allow for quick, versatile setups. When combined with S-Lock Quick Angle handles, the system gives surgeons control over each blade's exact placement while remaining extremely low profile. The system contains a streamlined selection of blades ideal for trauma procedures.




"When someone asks for an assistant, I say, 'Pull out a Thompson!'"

- Robert Mikkelsen, MD



Quick + Easy Blade Attachment

Add / remove blades, or switch from swivel to locked positions with the simple push of a button.




Quick Angle Technology

Quickly and precisely angle blades 45° up or down, while the low profile design keeps handles out of the way.



Versatile Frame Set Ups

Frame arms are independently adjusted, which helps achieve low profile retraction for a variety of procedures.



TRAUMA SYSTEM

#SL82036

REF	QTY	ITEM DESCRIPTION	PART #
A	1	Elite III Rail Clamp w/1 Cam Joint 18"	43902AC
B	1	Crossbar w/2 Cam III Joints 24" (11" x 13")	43900BC
C	2	Angled Arm 20" (10" x 10") @ 45°	44120
D	2	Angled Arm 24" (8" x 16") @ 45°	44124N
E	6	Cam II Clip-On Quick Angle 8"	SL42126WGP
F	2	T-Handle	60020
	1	Instrument Case 23" x 11" x 3 1/2"	50000G
G	2	Malleable 51mm x 165mm (2" x 6 1/2")	SL46112
	1	Malleable 38mm x 152mm (1 1/2" x 6")	SL46120B
H	1	Malleable 76mm x 254mm (3" x 10")	SO46124A

REF	QTY	ITEM DESCRIPTION	PART #
I	2	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
J	2	Balfour [sm lips] 70mm x 73mm (2 3/4" x 3")	SO46140C
K	2	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	SL46143
L	1	Richardson 51mm x 102mm (2" x 4")	SL46150
M	1	Richardson 51mm x 127mm (2" x 5")	SL46151
N	1	Richardson 51mm x 152mm (2" x 6")	SL46152
O	1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
P	1	St. Marks [lip] 64mm x 178mm (2 1/2" x 7")	SO46188
	1	Instrument Case 23" x 11" x 3 1/2"	50000G



Visit us online: thompsonsurgical.com

For a Free Trial Call **1-800-227-7543** *

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2020 Thompson Surgical Instruments, Inc.

® S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc. Patents: tpat.com

Rev A
052820
tmss

