


## GYN Oncology System SUGGESTED SET UP

### Bilateral GYN Oncology



Additional GYN set ups shown on page 40.

### GYN ONCOLOGY SYSTEM #SL82011

REF	PART #	ITEM DESCRIPTION	QTY
<b>Tray 1</b>			
A	43903AC	Elite III Rail Clamp with 2 Cam Joints 18"	2
B	41900F	Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4"	1
C	44112SP	Lower Abdominal Bar 14" x 6"	1
D	43900BC	24" Crossbar with 2 Cam Joints	1
	50000BL	Instrument Case (26" x 10" x 5")	1
<b>Tray 2</b>			
E	41900PH	34" Bilateral Crossbar Hinged (11 1/2" x 8 1/2" x 11 1/2")	1
F	43112C	Cam III Serrated Slide-On Joint	1
G	44118N	18" Angled Arm (8" x 10") @ 45°	2
H	44124W	24" Angled Arm (11" x 13") @ 45°	2
I	SL42126WGP	Cam II Clip-on Quick Angle 8"	6
J	SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2
	SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	2
K	60020	T-Handle	2
	50000G	Instrument Case (23" x 11" x 3 1/2")	1
<b>Tray 3</b>			
L	SL46115	Malleable 13mm x 64mm (1/2" x 2 1/2")	1
	SL46116	Malleable 13mm x 152mm (1/2" x 6")	1
	SL46117	Malleable 13mm x 203mm (1/2" x 8")	2
M	SL46118	Malleable 25mm x 203mm (1" x 8")	1

( GYN ONCOLOGY SYSTEM #SL82011 CONTINUED )

REF	PART #	ITEM DESCRIPTION	QTY
	SL46118A	Malleable 25mm x 254mm (1" x 10")	1
	SL46118P	Soft Malleable 25mm x 102mm (1" x 4")	1
	SL46122	Malleable 51mm x 203mm (2" x 8")	1
	SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	2
	SO46124A	Malleable 76mm x 254mm (3" x 10")	1
N	SO46124B	Malleable 102mm x 254mm (4" x 10")	1
	SL46127	Kelly 38mm x 51mm (1 1/2" x 2")	1
	SL46128	Kelly 51mm x 64mm (2" x 2 1/2")	1
	SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	1
O	SO46130	Kelly 76mm x 89mm (3 1/2" x 3")	2
P	SO46135	Fence 102mm x 127mm (4" x 5")	1
Q	SO46140C	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	1
	SO46141C	Balfour [small lips] 78mm x 73mm (3" x 3")	2
	SO46141D	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	2
	SL46142	Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	2
R	SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2
	SL46150	Richardson 51mm x 102mm (2" x 4")	1
S	SO46151	Richardson 51mm x 127mm (2" x 5")	1
	SO46153	Richardson 57mm x 178mm (2" x 7")	1
T	SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1
U	SO46170B	Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2")	1
V	SO46171B	Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2")	1
W	SL46174	Deaver 51mm x 127mm (2" x 5")	1
X	SL46180	Heany 25mm x 102mm (1" x 4")	1
Y	SO46188	St. Mark's 64mm x 178mm (2 1/2" x 7")	1
Z	SL46191	Malleable Finger 152mm (6")	2
	SO46241A	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	1
1	SO46241B	Wide Balfour 100mm x 61mm (4" x 2 3/8")	1
	SO46241C	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	2
2	SL46660	Double Angled 25mm x 76mm (1" x 3")	2
3	SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2
	50000G	Instrument Case (23" x 11" x 3 1/2")	1

# GYN Oncology System COMPONENTS

RETRACTOR SYSTEMS

