


**NEW Trauma System** SUGGESTED SET UP

**Abdominal Trauma**



“When someone asks for an assistant,  
I say, ‘Pull out a Thompson!’”

- Robert Mikkelsen, MD

**TRAUMA SYSTEM**


**#SL82036**

| REF           | PART #     | ITEM DESCRIPTION                               | QTY |
|---------------|------------|--|-----|
| <b>Tray 1</b> |            |  |     |
| A             | 43902AC    | Elite III Rail Clamp with 1 Cam Joint 18"      | 1   |
| B             | 43900BC    | 24" Crossbar with 2 Cam III Joints (11" x 13") | 1   |
| C             | 44120      | 20" Angled Arm (10" x 10") @ 45°               | 2   |
| D             | 44124N     | 24" Angled Arm (8" x 16") @ 45°                | 2   |
| E             | SL42126WGP | Cam II Clip-On Quick Angle 8"                  | 6   |
| F             | 60020      | T-Handle                                       | 2   |
|               | 50000G     | Instrument Case (23" x 11" x 3 1/2")           | 1   |
| <b>Tray 2</b> |            |  |     |
| G             | SL46112    | Malleable 51mm x 165mm (2" x 6 1/2")           | 2   |
|               | SL46120B   | Malleable 38mm x 152mm (1 1/2" x 6")           | 1   |
| H             | SO46124A   | Malleable 76mm x 254mm (3" x 10")              | 1   |
| I             | SO46129    | Kelly 64mm x 76mm (2 1/2" x 3")                | 2   |
| J             | SO46140C   | Balfour [sm lips] 70mm x 73mm (2 3/4" x 3")    | 2   |
| K             | SL46143    | Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")      | 2   |
| L             | SL46150    | Richardson 51mm x 102mm (2" x 4")              | 1   |
| M             | SL46151    | Richardson 51mm x 127mm (2" x 5")              | 1   |
| N             | SL46152    | Richardson 51mm x 152mm (2" x 6")              | 1   |
| O             | SL46160    | Harrington 64mm x 152mm (2 1/2" x 6")          | 1   |
| P             | SO46188    | St. Marks [lip] 64mm x 178mm (2 1/2" x 7")     | 1   |
|               | 50000G     | Instrument Case (23" x 11" x 3 1/2")           | 1   |

### Trauma System COMPONENTS



**NEW** Trauma System contains a streamlined selection of blades ideal for trauma procedures.



RETRACTOR SYSTEMS