

"This system is a shoulder-saver and is like having another tireless assistant standing next to me."

- Steven Thomas, MD

Breast Retractor System

The Thompson Breast Retractor System allows for a quick set-up, virtually unlimited positioning, and rigid, stable retraction where you need it most. The large selection of blades accommodates breast procedures such as mastectomies and breast reconstruction, all while suiting varying patient anatomies.




S-Lock® Technology

Swivel 360° or lock in any position. Switch from swivel to locked with the push of a button.



Fast, Easily Adjustable Arms


Multi-jointed Articulating Arms allow for easy repositioning of blades and multi-planar retraction.



Versatile Set Up Options

S-Lock Articulating Arms offer simple, quick retraction, while the C-Frame utilizes handles with additional features.

Breast System ORDERING INFORMATION



BREAST SYSTEM

#SL82041

REF	PART #	ITEM DESCRIPTION	QTY
Tray 1			
A	43903AC	Elite III Rail Clamp w/2 Cam Joints 18"	1
B	SL42191	S-Lock Articulating Arm Small 15 1/2"	2
C	SL46126	Richardson 25mm x 32mm (1" x 1 1/4")	2
D	SL46422ET	Radio, Richardson 50mm x 60mm (2" x 2 3/8")	2
E	SL46426ET	Radio, Richardson 50mm x 100mm (2" x 4")	2
F	SL46432ET	Radio, Richardson 50mm x 160mm (2" x 6 1/4")	2
	SL46437ET	Radiolucent Richardson 50mm x 210mm (2" x 8 1/4")	2
	SL46441ET	Radiolucent Richardson 50mm x 250mm (2" x 9 7/8")	2
G	SO46402	Cat's Paw 6mm x 38mm (1/4" x 1 1/2")	2
H	SL46214	Murphy Sharp 4 Prong 25mm x 13mm (1" x 1/2")	2
I	SL46115	Malleable 13mm x 64mm (1/2" x 2 1/2")	1

REF	PART #	ITEM DESCRIPTION	QTY
J	SL46116	Malleable 13mm x 152mm (1/2" x 6")	1
K	SL46119	Malleable 25mm x 127mm (1" x 5")	1
	50000G	Instrument Case (23"x 11" x 3 1/2")	1
Tray 2			
L	44615FH	14" Extension Arm with Fixed Hinge Joint	1
M	44225	C-Frame 13 1/2"	1
N	SL42126WGLP	Cam II Clip-on Quick Angle 10"	3
O	SL45011CGP	Mini Micro-Adjustable II Clip-on QA 9"	3
P	60020	T-Handle	2
	50000G	Instrument Case (23"x 11" x 3 1/2")	1
Tray 3			
Q	40020X	Lite Wand Xe	1
	50000LW	Instrument Tray - Lite Wand	1